

# NanoHub

# NanoHub

### The smallest touchscreen anesthesia Tabletop station

Delivered with one induction chamber and up to two masks. Dedicated for rodents, with possible application on birds and all other small animals





- User safety protection (less than 2 ppm),
- Intuitive touchscreen use,
- Accurate independant gas flow for one chamber and two masks,
- Isoflurane and sevoflurane compatible.
- Easy connection to carrier gas (air and/or oxygen)

### Technical specifications.

#### Gas management

Electronic regulation and touchscreen software operating. TEC3 vaporizer with keyFillTM secured filling system. Adjustable flow for rodents,

Inlet/outlet for one induction and two independant masks.

#### Induction

One induction chamber with electronic induction Timer. Opening/closing safety system.



One or two masks with large choice of animal species (rat / mice / rabbit ).

#### Scavenging system

Pneumatic scavenger with two independent channels. One channel for induction chambers, one channel for masks. F/AIR canisters to filter halogen used gas. Filter saturation indicator display.

#### TEMSEGA : your designer-manufacturer for all-in-one animal anesthesia solutions.

Since 1992, TEM SEGA is the leading manufacturer in Europe of gaseous anesthesia devices for veterinarian purpose. Our device can be customized to fit many species from mice to horses, and drive all types of labs gas (air, oxygen, CO2, Nitrogen and toxic gas). Our values are to offer the highest human and animal protection, improve productivity and precision in research labs, propose modularity and flexibility, and comply with ethic guidelines. Our technical team install, train and support everywhere in the world.

\*Our equipement are complying with European and US regulation, AAALAC recommendations.







## **Accessories & options**

Other vaporizers available: TEC850, second-hand TEC3. Stereotaxic mask, induction chamber for each model of imager.



## Informations

Find informations about our products on our website : https://temsega.com



27 allée Charles Darwin 33600 PESSAC +33 5 56 32 29 63 info@temsega.fr www.temsega.com